

**GLANS COMPATIBLE SINGLE UNIT SEMEN COLLECTION AND  
STORAGE DEVICE, KIT AND RELATED METHODS OF USE**

**5 ABSTRACT OF THE DISCLOSURE**

A method and device to conform with the glans penis to recover ejaculated semen completely, to prevent the loss of initial sperm rich epididymal fractions, to avoid the use of a condom for masturbation, to eliminate the multi-step transfers of semen following ejaculation that are common with current methods, and to provide a  
10 single device that contains a semen collecting chamber that fits onto a glans penis. The device provides a storage and measuring reservoir, an optional cap, and a vertical stand in one integrally formed module in a sterile pack. The device and method will have multiple uses in a variety of contexts including, but not limited thereto: in the diagnosis of infertility, in semen donation, in vitro fertilization (IVF)/ intra-  
15 cytoplasmic sperm injection (ICSI) clinics, hospitals, forensic laboratories and research laboratories and will be included in kits intended for over-the-counter sperm testing devices

20